



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Rubber Cork for Hypsometer

Product Code :

SCL-304/050

Image not found or type unknown



Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Rubber Cork for Hypsometer (Product Code: SCL-304/050) is a red rubber bung used to seal the thermometer and tube openings of a hypsometer during thermometry experiments in physics laboratories. Lab staff order it as a spare when the original bungs harden, crack or go missing, and it is offered as an assorted lot. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, from its range of physics instruments for schools and colleges.

Key Details

Parameter	Detail
Item	Rubber cork (bung) for hypsometer
Material	Rubber
Colour	Red
Assortment	Assorted; the mix is confirmed with each quotation
Product code	SCL-304/050
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is the Rubber Cork for Hypsometer Used For?

A hypsometer depends on tight bungs: one holds the thermometer under test in the lid tubule and another seals the manometric tube in its side tubulure. If a bung is cracked or loose, steam escapes and both the pressure check and the fixed-point reading lose their meaning. These rubber corks replace worn or lost bungs so the apparatus can go straight back into use.

When a hole is needed, bore it with a sharp cork borer of the right size. Wet the thermometer with a little water or glycerine, hold it close to the cork with a cloth and ease it through with a gentle twist; never force glass through a tight hole. Keep spare corks away from heat and oil, which make rubber harden or swell.

What to Specify When Ordering

- The openings the corks must fit: the lid tubule for the thermometer, the manometric tube neck, or both
- The number of corks needed per hypsometer and per lab
- Whether a complete Hypsometer or other physics apparatus should be quoted in the same order

Supplier

Jain Scientific Equipments Private Limited trades as ScienceLab and supplies these rubber corks from Ambala, Haryana, as spares for the physics apparatus it offers schools and colleges. Ordering corks together with the hypsometer keeps a thermometry set complete, and the same enquiry can cover other consumables. Company details are on the [ScienceLab company profile](#).

Related

Useful co

Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803