



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Geiger Muller Tube Holder

Product Code :

SLE/PHY/42



Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Geiger Muller Tube Holder (Product Code: SLE/PHY/42) is an aluminium box housing a G.M. tube socket, used to mount a Geiger–Müller tube and connect it to a counter for radioactivity measurements in physics laboratories. Designed for use with the Scaler Timer, it connects through a 1 m co-axial lead ending in a standard P.E.T. connector and has a stabilising resistor already wired inside. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India.

Specifications

Parameter	Detail
Housing	Aluminium box containing the G.M. tube socket
Connecting lead	1 m long co-axial lead
Connector	Standard P.E.T. connector
Anode resistor	2.2 megohm stabilising resistor, wired internally in series with the anode connection
For use with	Scaler Timer
Product code	SLE/PHY/42
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is the GM Tube Holder Used For?

The Geiger Muller Tube Holder turns a bare GM tube into a bench-ready detector. The tube plugs into the socket inside the aluminium box, and the co-axial lead links it to the counter, so students can place the detector at a measured distance from a source and record counts. Common practicals include background counts, count rate against distance, absorption by paper or metal sheets, and half-life measurement.

Because the 2.2 megohm stabilising resistor sits inside the box in series with the anode, there is no loose resistor to wire on the bench, which keeps a Geiger Muller tube set-up in an Indian school or college lab tidy and repeatable. Before counting, set the tube voltage on the counter to suit your GM tube, and keep the delicate tube window away from fingers and sharp objects.

What to Specify When Ordering

- Number of detector stations, i.e. how many Geiger Muller Tube Holders you need.
- Whether you also need the Scaler Timer the holder is designed to work with.
- Whether you need a GM tube, as a tube is not listed as part of the holder.

Supplier

ScienceLab, the trading name of Jain Scientific Equipments Private Limited, supplies the Geiger Muller Tube Holder along with counting instruments and radioactivity accessories for physics laboratories. For more information, please contact us at sales@scilabexport.com, quoted together with your requirements.

Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

Related

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803