



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Laboratory digital planetary mixer

Product Code :

SLE/CE/75



Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Laboratory digital planetary mixer (Product Code: SLE/CE/75) is a table-mounted planetary mixer used to blend asphalt mixes evenly for specimen preparation in highway and materials testing laboratories. Road research labs and engineering colleges use its 5 litre bowl and whisk with user-set planetary speeds of 5 to 65 rpm, and a front grill that stops the machine as soon as it is opened. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, a laboratory equipment supplier for highway and materials labs.

Specifications

Parameter	Detail
Application	Mixing asphalt mixes; can be used as part of a Superpave performance-based testing programme
Mounting and action	Table-mounted unit with planetary mixing action
Bowl capacity	5 litres
Supplied with	Bowl and whisk, both easily fitted and removed
Planetary speed	User-defined, 5 to 65 rpm
Whisk speed	User-defined, 12 to 149 rpm
Speed control	Easily selected and adjustable during mixing
Controls	Dedicated display and keyboard interface; alphanumeric display, 2 x 16 characters
Safety	Opening the front grill automatically stops the machine
Overall dimensions	465 x 540 x 620 mm (l x d x h)
Weight	Approx. 35 kg
Power	370 W
Product code	SLE/CE/75
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is the Digital Planetary Mixer Used For?

Every trial blend in an asphalt mix design has to be mixed thoroughly so the binder coats all the aggregate before specimens are compacted, and doing that by hand in a pan is slow and varies from person to person. In a planetary mixer the whisk spins on its own axis while also travelling around the bowl, so it reaches the whole batch. The operator keys in planetary and whisk speeds on the display and can change them while the batch is running.

Mixed batches then move on to compaction and to tests such as stability and maximum density, including within a Superpave performance-based programme. For buyers comparing laboratory instruments, the digital planetary mixer will be a time-saving investment for materials labs. Weigh balance is also a common requirement for materials testing.

Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

What to 1226/1-5, Bengali Mohalla, Science Market | 133001, Hariana



Email : sales@scilabexport.com

Phone: +91-7082934803