



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Magnetic Stirrer Bar

Product Code :

SLE/LS/10



Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Magnetic Stirrer Bar (Product Code: SLE/LS/10) is a small magnetic bar placed inside a beaker or flask so that a magnetic stirrer can spin it and mix the liquid, used in chemistry, biology and research laboratories. Each bar measures 25 × 6 mm, and they come in a pack of 10, giving every stirrer or student group a spare. They are supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, as lab equipment accessories for buyers in India and abroad.

Specifications

Parameter	Detail
Dimensions	25 × 6 mm
Pack size	10 bars
Use	Stirring element (follower) for magnetic stirrers and hotplate stirrers
Product code	SLE/LS/10
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is a Magnetic Stirrer Bar Used For?

The bar — also called a follower or flea — follows the rotating magnetic field under the stirrerplate and turns inside the vessel, mixing a solution without a glass rod or an overhead stirrer. Bars are easily lost down sinks, left in waste flasks or scratched with use, so a pack of spares keeps practicals running. Departments often order them with other consumables from one lab consumables supplier in India.

Before ordering, check that your stirrer accepts a bar of this length and that the bar lies flat across the base of your beakers and flasks. Take bars out with a retriever rather than fingers, rinse and dry them, and set aside any bar whose surface is cracked or chipped.

What to Specify When Ordering

- The number of packs of 10 needed
- That the 25 × 6 mm size suits your stirrers and vessels
- Whether a stirrer bar retriever should be added to the order

Supplier

Jain Scientific Equipments Private Limited trades as ScienceLab and supplies these stirrer bars from Ambala, Haryana, alongside the magnetic stirrers and hotplate stirrers in its range. Adding packs to a stirrer order means the spares arrive with the instruments. Learn more on our [company profile](#).

Related Laboratory Stirrer

Stirrers and tools used with these bars:

- [Magnetic Stirrer](#)
- [Magnetic Stirrer](#)
- [Magnetic Stirrer](#)
- [Stirrer](#)

Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803