



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Mariot Bottle (for distilled water)

Product Code :

SCI-612/102

Image not found or type unknown



Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Mariot Bottle (for distilled water) (Product Code: SCI-612/102) is a 2.5 litre transparent glass bottle with a tap, used to store distilled water on the bench and dispense it as needed in school, college and routine laboratories. The tap lets students draw water into beakers or wash bottles without lifting or tipping the full bottle. It is manufactured by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, an OEM manufacturer, at its own factory in Ambala, India.

Specifications

Parameter	Detail
Capacity	2.5 litre
Outlet	With tap
Material	Glass
Colour	Transparent
Chemical resistance	Chemical resistant
Product code	SCI-612/102
Manufacturer	Jain Scientific Equipments Pvt Ltd (ScienceLab) — OEM
Country of origin	India (made in Ambala, Haryana)

What Is the Mariot (Mariotte) Bottle Used For?

Practicals that need a steady supply of distilled water, such as making up solutions, rinsing glassware or filling burettes, run more smoothly with a tapped reservoir at the end of the bench. The 2.5 litre bottle is filled from the still or store container and set on a raised shelf or stand, and students draw water through the tap into a beaker or wash bottle.

The name refers to the Mariotte bottle, a classic laboratory design for delivering liquid at a steady rate. If your experiment needs that constant-flow arrangement, ask us to confirm the stopper and tube details of the bottle you are quoted. Keep the bottle closed to limit dust and carbon dioxide pick-up from the air, and rinse it with fresh distilled water at each refill.

What to Specify When Ordering

- Bottle: 2.5 litre with tap, code SCI-612/102.
- Flow arrangement: whether you need a constant-flow (Mariotte) stopper and tube set-up — confirm with us.
- Placement: a raised shelf or stand so the tap clears beakers and wash bottles.
- Quantity: typically one reservoir per lab or bench block.

OEM Manufacturer

Jain Scientific Equipments Private Limited (CIN U29309HR2020PTC087597) trades as ScienceLab and is the OEM of the Mariot Bottle (for distilled water). The bottle is made at the company's factory in Ambala, Haryana.

Licence A ScienceLab is available for sale to schools, colleges, and laboratories. For more information, please contact us. We are also available for OEM orders. **Website:** <https://www.scilabexport.com>, **Email:** sales@scilabexport.com
1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803