



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803


**Product Name :**

Brush

**Product Code :**

"SCL282/046"

Image not found or type unknown



**Website:** <https://www.scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

## Description :

**Brush** (Product Code: SCL282/046) is a plastic-bodied cleaning brush used to scrub glassware and apparatus clean after practicals in school, college and training laboratories. Lab assistants keep it by the sink because food tests, stains and culture work leave films that a rinse alone will not shift. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, as part of its range as a lab consumables supplier in India.

## Key Details

| Parameter     | Detail                                                         |
|---------------|----------------------------------------------------------------|
| Type          | Laboratory cleaning brush                                      |
| Body material | Plastic                                                        |
| Purpose       | Used for cleaning                                              |
| Product code  | SCL282/046                                                     |
| Supplied by   | Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India |

## What Is the Lab Cleaning Brush Used For?

After a session, empty the apparatus, rinse out loose material, and then work warm water and laboratory detergent over the inside surfaces with the brush. Pay attention to the base and shoulders of vessels, where precipitates from Benedict's tests or dried stain collect. Finish with a tap-water rinse and a final rinse in distilled water, then leave items upside down to drain.

Good cleaning habits are an everyday part of lab equipment maintenance in India's school laboratories and prevent carry-over that can spoil the next group's results. Use gentle pressure on thin glassware and never force the brush into a neck that is too narrow for it. It also helps to keep one brush for chemical residues and another for biological material, labelled so they are not mixed up.

## Supplier

ScienceLab is the trading name of **Jain Scientific Equipments Private Limited**, an Ambala, Haryana company that supplies this brush together with glassware, apparatus and consumables for educational laboratories. Ordering cleaning supplies with the rest of a lab list means one quotation and one delivery for the whole practical set-up. See the [company profile](#) for more.

## Related Laboratory Equipment

Apparatus this brush typically cleans, plus shaped brushes and rinse water:

- [Test Tubes \(With Rim\)](#) — the items most often scrubbed after food-test practicals.
- [Beaker](#) — collects stains and precipitates that need brushing out.
- [Brush Test Tube](#) — a tube-shaped brush for narrow tubes.
- [Beaker Brush](#) — a brush suited to wide-mouthed beakers.
- [Distilled Water](#) — for the final rinse so no tap water residue remains.

The full B

**Website:** <https://www.scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com)

Frequen

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803