



Email : sales@scilabexport.com

Phone: +91-7082934803


Product Name :

Test Tube Brush Nylon

Product Code :

SCL-630/063

Image not found or type unknown



Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market, Ambala Cantt. Harvana 133001, India



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

The **Test Tube Brush Nylon** (Product Code: SCL-630/063) is a round cleaning brush with soft, white nylon bristles, made for scrubbing precipitates, stains and dried residues from the inside of test tubes after practical work. The soft bristles clean the curved bottom and walls of the tube without the harshness of abrasive scouring. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, a supplier of chemistry consumables and chemistry lab equipment in India to schools, colleges and laboratory stores.

Specifications

Parameter	Detail
Application	Cleaning test tubes
Bristle material	Nylon
Bristle stiffness	Soft
Shape	Round
Colour	White

Typical Use

Clean glassware is the starting point of reliable results — a trace of old precipitate can spoil a salt-analysis test, and the silver-mirror test in organic chemistry only works in a thoroughly clean tube. At the end of each practical, empty the tube, rinse it, add a little detergent solution and work the brush up and down with a twisting motion so the bristles reach the rounded bottom. Rinse with tap water, give a final rinse with distilled water from a wash bottle, and stand the tube upside down to drain.

Push the brush in gently rather than jabbing it against the base of the tube, which protects thin glass from knocks. Rinse and shake the brush after use and hang it to dry so the nylon stays fresh; replace it once the bristles become flattened or splayed.

Supplier

This brush is supplied by **Jain Scientific Equipments Private Limited**, which trades as **ScienceLab** from Ambala, Haryana, India. For vendor registration, the company's details are GSTIN 06AAECJ8618A1ZT and Importer-Exporter Code AAECJ8618A issued by the DGFT, under which it ships export orders. Find out more on the [company profile](#).

Related Laboratory Equipment

Other items for the glassware washing routine:

- [Test Tube Stand Poly Lab](#) — to hold cleaned tubes while they drain
- [Boiling Tubes](#) — heavier tubes that need the same regular cleaning
- [Brush Burette](#) — a longer brush for cleaning burettes
- [Wash Bottle Burette](#) — a longer brush for cleaning burettes
- [Wash Bottle Burette](#) — a longer brush for cleaning burettes

More

Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market, Ambala Cantt. Harvana 133001, India



Email : sales@scilabexport.com

Phone: +91-7082934803

Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market, Ambala Cantt. Harvana 133001, India