



Email : sales@scilabexport.com

Phone: +91-7082934803


Product Name :

Methylated Spirit

Product Code :

SCL-896/191

Image not found or type unknown



Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market, Ambala Cantt. Harvana 133001, India



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Methylated Spirit (Product Code: SCL-896/191) is denatured ethyl alcohol used as spirit lamp fuel and as a cleaning and rinsing solvent in school, college and teaching laboratories. It is packed in a 5 litre stock container of analytical reagent grade with an MSDS, suiting labs that heat and clean glassware every working day. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, a laboratory chemicals supplier in India, as a standing item on the lab consumables list of institutions across the country.

Specifications

Parameter	Detail
Product	Methylated spirit (denatured ethyl alcohol)
Capacity (pack size)	5 litre
Grade	Analytical reagent
Safety documentation	Supplied with MSDS (Material Safety Data Sheet)
Shelf life	Shelf life depends on the chemical, its grade and storage conditions; ask us for the expiry date of the batch you are quoted.
Product code	SCL-896/191
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Is Methylated Spirit Used For?

In the school laboratory, methylated spirit is the usual fuel for glass spirit lamps, giving a steady flame for gentle heating of test tubes and for flaming inoculating loops where no gas line is available. It is also used to wipe down benches, to rinse the last traces of water from freshly washed glassware so that it dries quickly, and to dissolve grease, resin and marker ink that water alone will not remove.

The 5 litre container is a stock size, so decant a working quantity into a labelled wash bottle and refill lamps from that smaller bottle instead of pouring from the large can near students. Fill lamps only when they are cold and unlit, keep the container closed and away from flames and sunlight, and read the MSDS before first use. Schools planning their school chemistry lab equipment in India usually keep methylated spirit on the standing consumables list next to burners and spare wicks.


Supplier

Jain Scientific Equipments Private Limited trades under the name **ScienceLab** and is based in Ambala, Haryana, India. The company supplies chemicals such as this methylated spirit to schools, colleges and dealers in India and international markets, so that an institution can source apparatus and consumables from one place. You can read more on our [company profile](#) page.

Related Chemistry Lab Products

Items **Website:** <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

- 1226/1-5, Bengali Mohalla, Science Market, Ambala Cantt, Harvana 133001, India



Email : sales@scilabexport.com

Phone: +91-7082934803

Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market, Ambala Cantt. Harvana 133001, India