



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Beaker Brush

Product Code :

SLE/GB/11



Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Beaker Brush (Product Code: SLE/GB/11) is a cleaning brush for beakers, with a 4 inch section of nylon bristles on an iron body, 14 inch long overall and 20 mm in diameter. The long stem reaches the base of large and tall beakers, while the nylon bristles scrub away dried residues, precipitates and indicator stains from the walls. It is supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, as one of the everyday items it offers as a lab consumables supplier in India.

Specifications

Parameter	Detail
Purpose	Cleaning beakers
Bristle material	Nylon
Bristle length	4 inch
Body material	Iron
Overall length	14 inch
Diameter	20 mm

Typical Use

Clean glassware is the starting point of reliable results, so a beaker brush sits at every lab sink. After a practical, rinse out the beaker, add warm water with a little laboratory detergent, scrub the base and walls with the brush using a twisting motion, then rinse with tap water followed by distilled water and leave the beaker upside down on a drying rack. The 4 inch bristle head covers a large area per stroke, which speeds up washing after a full class has used the glassware.

Avoid pushing the bare end of the stem against the glass, as metal can scratch or chip it. Rinse the brush after use and let it dry before storing, which keeps the iron body from rusting and the bristles in shape; replace it once the bristles splay or wear down.

Supplier

This brush is supplied by Jain Scientific Equipments Private Limited, which trades as ScienceLab from Ambala, Haryana, as part of the cleaning and consumable items it offers alongside its laboratory apparatus. For vendor registration, the company's details are GSTIN 06AAECJ8618A1ZT and Importer-Exporter Code AAECJ8618A from the DGFT, under which its export orders are handled. More is available in our [company profile](#).

Related Laboratory Equipment

Other items for a well-equipped glassware washing station:

- [Brush, Beaker, Nylon](#) — an alternative nylon beaker brush
- [Bottle Brush](#) — for reagent and narrow neck bottles
- [Brush, Beaker, Nylon](#)
- [Washing Station](#)
- [Drying Rack](#)

Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803