



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

**Product Name :**

Widemouthed Wash Bottles

**Product Code :**

SLE/P/03





Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803

## Description :

**Widemouthed Wash Bottles** (Product Code: SLE/P/03) are 500 mL (16 oz) squeeze bottles used to rinse glassware, wash down precipitates and add liquid drop by drop in chemistry, biology and clinical labs. Students and technicians like them because the wide mouth makes refilling quick and clean, and the standard spout delivers either a single drop or a needle-sharp stream. They are supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, laboratory equipment supplier in India for teaching and routine labs.

## Specifications

| Parameter            | Detail                                                                       |
|----------------------|------------------------------------------------------------------------------|
| Capacity             | 500 mL (16 oz)                                                               |
| Mouth                | Wide, for easy filling                                                       |
| Spout style          | Standard spout                                                               |
| Delivery             | Single drop or needle-sharp stream; tip can be cut back to increase delivery |
| Material (as listed) | LDPE / polypropylene                                                         |
| Product code         | SLE/P/03                                                                     |
| Supplied by          | Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India               |

## What Are Widemouthed Wash Bottles Used For?

A wash bottle is the most-used tool on a chemistry bench. It rinses the last traces of a solution from a beaker into a volumetric flask, washes a precipitate down the side of a filter funnel, and brings a flask up to the mark one drop at a time. The needle-sharp stream reaches corners of glassware that a tap cannot, while a gentle squeeze gives a single controlled drop.

The wide mouth is the practical advantage of this bottle: it can be refilled straight from a carboy or jug without a funnel, and its inside is easy to clean when the contents change. Label each bottle with what it holds, such as distilled water or ethanol, and cut the tip back only a little at a time, because a trimmed tip cannot be made finer again.

## Supplier

Jain Scientific Equipments Private Limited supplies the Widemouthed Wash Bottles from Ambala, Haryana, under its ScienceLab trading name, as part of the Scientific Lab Plasticware range. Schools, colleges, pathology and industrial labs, dealers and importers in India and international markets can order them with beakers, bottles and funnels on one enquiry. More about the company is on the [profile page](#).

## Related Scientific Lab Plasticware

Other rinsing and filtration items that work with this wash bottle:

- [Un](#)
  - [Po](#)
  - [Wa](#)
- Website:** <https://www.scilabexport.com>, **Email:** [sales@scilabexport.com](mailto:sales@scilabexport.com) one size.
- 1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana ter a



Email : [sales@scilabexport.com](mailto:sales@scilabexport.com)

Phone: +91-7082934803