



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Unitary Wash Bottles

Product Code :

SLE/P/05



Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Unitary Wash Bottles (Product Code: SLE/P/05) are plastic wash bottles whose bottle and delivery tubing are moulded in one piece, used to rinse glassware and dispense wash liquids in teaching, research and clinical labs. Because there is no joint between bottle and tube, there is nothing to work loose or leak, and the whole contents can be dispensed without tipping or shaking. They are supplied by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, Ambala, India, a school lab equipment supplier in India, in three listed sizes.

Specifications

Parameter	Detail
Construction	Bottle and tubing moulded in one piece for strength, with no joint to leak
Dispensing	Entire contents can be dispensed without tipping or shaking
Tip	Delivers a fine stream; can be cut back to increase flow
Material	Plastic
Sizes (as listed)	125, 500 and 1000 (unit not stated in the listing; confirm when ordering)
Product code	SLE/P/05
Supplied by	Jain Scientific Equipments Pvt Ltd (ScienceLab), Ambala, India

What Are Unitary Wash Bottles Used For?

On a busy student bench, the weak point of an ordinary wash bottle is the joint where the tube passes through the cap: it loosens, drips and lets liquid run down the hand. The one-piece design of the Unitary Wash Bottles removes that joint, so the bottle can be squeezed hard, knocked over or passed around a class without leaking. And since the whole contents can be dispensed without tipping, the bottle stays upright until the last of the liquid is used.

Use the fine stream for rinsing burettes, pipettes and the walls of titration flasks, and cut the tip back a little if you need faster flow for washing large glassware. Keep separate, clearly labelled bottles for distilled water, ethanol or other wash liquids so they are never mixed up.

What to Specify When Ordering

- Which listed size you need — 125, 500 or 1000 — and confirm the unit with us.
- The number of bottles of each size.
- The liquids you plan to use, so we can check they suit the plastic.

Supplier

ScienceLab, the trading name of Jain Scientific Equipments Private Limited in Ambala, Haryana, supplies the Unitary Wash Bottles within its Scientific Lab Plasticware range. They can be ordered by schools, universities, hospital and industrial labs, dealers and importers in India and international markets, on their own or with other bench plastic ware. See the [company profile](#) for more.

Related **Website:** <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803