



Email : sales@scilabexport.com


Phone: +91-7082934803

Product Name :

Kamje Pirosti

Product Code :

SCL324/062



Science Lab



Email : sales@scilabexport.com
Phone: +91-7082934803

Description :

Kamje Pirosti (Product Code: SCL324/062) is a triangular laboratory tripod stand, approximately 200 mm high, used to support a gauze and a beaker or flask above a spirit lamp or burner during heating in school laboratories. Its three-legged triangular frame gives a steady base on the bench, making it a core piece whenever a tripod stand for a lab in India is on the list. It is made by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, an OEM manufacturer, at its own factory in Ambala, India.

Key Details

Parameter	Detail
Product type	Laboratory tripod stand
Shape	Triangular
Height	Approximately 200 mm
Product code	SCL324/062
Manufacturer	Jain Scientific Equipments Pvt Ltd (ScienceLab) — OEM
Country of origin	India (made in Ambala, Haryana)

What Is the Tripod Stand Used For?

The tripod is the frame of the classic school heating set-up: a gauze rests across its triangular frame, a beaker or conical flask stands on the gauze, and a spirit lamp or Bunsen burner goes underneath. Students use it to boil water, dissolve solids faster, evaporate a salt solution or warm a water bath, keeping the vessel steady and at a safe height above the flame.

Before the practical, check that the tip of your lamp or burner flame reaches just below the gauze at this roughly 200 mm height. Stand the tripod on a heat-resistant mat, keep it clear of the bench edge, and remember that the stand stays hot for some time after the flame is out, so move it only once it has cooled.

OEM Manufacturer

Kamje Pirosti is manufactured by ScienceLab, which is the OEM for this stand. Jain Scientific Equipments Private Limited (CIN U29309HR2020PTC087597), the company behind the ScienceLab brand, carries out production at its licensed factory in Saha, Ambala (Factory Licence AMB-ONLINE-CHD-J-414). Schools, colleges and distributors therefore buy directly from the people who make the stand, and quotations and repeat orders are dealt with by the manufacturer itself. Our [company profile](#) has more information.

Related Laboratory Equipment

The other parts of a bench heating station:

- [Ceramic Mesh](#) 125 x 125 mm, used to support the tripod
- [Alcohol Lamp](#)
- [Alcohol Burner](#)
- [Cotton Gauze](#)
- [Tri](#)

Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana



Email : sales@scilabexport.com

Phone: +91-7082934803