



Email : sales@scilabexport.com

Phone: +91-7082934803

Product Name :

Clock Glass 45mm, 60mm, 70mm

Product Code :

SCI-601/106

Image not found or type unknown



Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Harvana



Email : sales@scilabexport.com

Phone: +91-7082934803

Description :

Clock Glass 45mm, 60mm, 70mm (Product Code: SCI-601/106) is a shallow, curved disc of borosilicate glass, known in most laboratories as a watch glass, used to hold, weigh, cover or evaporate small quantities of material in chemistry and biology practicals. It comes in three diameters — 45 mm, 60 mm and 70 mm — for students and technicians alike, and its smooth glazed finish keeps residues easy to see. It is made by **Jain Scientific Equipments Pvt Ltd (ScienceLab)**, an OEM manufacturer, at its own factory in Ambala, India.

Specifications

Parameter	Detail
Diameters available	45 mm, 60 mm, 70 mm
Material	Borosilicate glass
Finish	Glazed, smooth
Clarity	Clear, for easy viewing of the contents
Build	Durable
Also known as	Watch glass
Product code	SCI-601/106
Manufacturer	Jain Scientific Equipments Pvt Ltd (ScienceLab) — OEM
Country of origin	India (made in Ambala, Haryana)

What Is the Clock Glass (Watch Glass) Used For?

A clock glass is one of the most versatile small items on a chemistry bench. Students place a few crystals on it to weigh them on a balance, let a few drops of a salt solution evaporate on it to watch crystals appear, or lay it over a beaker to keep dust out while a solution stands. In biology it holds a small specimen, a few seeds or a drop of stain for inspection with a hand lens.

Pick a diameter slightly larger than the rim of the beaker you want to cover, and keep the smallest size for spot tests with indicators. Borosilicate glass tolerates gentle warming, but warm the glass over a water bath or on a gauze rather than in a direct flame, and let it cool before washing. Schools building up a stock of borosilicate laboratory glassware in India often order the three sizes together.

What to Specify When Ordering

- Diameter: 45 mm, 60 mm or 70 mm — list each size you need.
- Quantity per diameter, for example one per student pair or per bench.
- The beakers you plan to cover, so the right diameter can be matched.
- Product code SCI-601/106 on your enquiry.

OEM Manufacturer

Jain Scientific Equipments Pvt Ltd is an OEM manufacturer of laboratory glassware. The

Website: <https://www.scilabexport.com>, **Email:** sales@scilabexport.com

1226/1-5, Bengali Mohalla, Science Market | 133001, Haryana

MSME
60mm,
and



Email : sales@scilabexport.com

Phone: +91-7082934803